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※For use with Labospect008 and 008α, please use QUALIGENT N Calibrator.

## Preparation

Add 2.0 mL of purified water to a vial of lyophilized calibrator. Cap the vial securely and allow it to stand for a few minutes at room temperature. In order to dissolve the lyophilized calibrator completely, gently invert the vial about 20 times. To measure CHO, TG, F-CHO, or PL, allow the vial to stand for at least 30 minutes. To measure HDL-C or LDL-C, allow the vial to stand for at least 1 hour. Immediately prior to use, gently invert the vial once again.

Please note that the assigned values for each test, printed on this kit insert, vary by lot.

The resulting calibrator solution is stable for 1 week after reconstitution when stored at 2-10°C.

## Handling Precautions

Human material used in formulating this product has been tested and found to be negative for HBs antigens, anti-HCV antibodies, and anti-HIV antibodies. However, handle with care, as transmission of infectious disease may still be possible. Please wear protective gloves when preparing this product.

## Packaging

3 vials of 2.0 mL (lyophilized)

## Storage Temperature

2-10°C

## Expiration Date

Refer to the external packaging.

Lot No. **1 3 7 R H U**

## Assigned value

Test	Product	Concentration (mg/dL)	expanded uncertainty (mg/dL)
HDL-cholesterol	Cholestest N HDL / QUALIGENT HDL	55. 0	± 1. 1
LDL-cholesterol	Cholestest LDL / QUALIGENT LDL	111. 0	± 2. 8
Total cholesterol	Pureauto S CHO-N	177	± 2. 3
	Cholestest CHO / QUALIGENT CHO	177	± 2. 3
Free cholesterol	Pureauto S F-CHO-N	42	—
Phospholipid	Pureauto S PL	200	—
Triglyceride	Pureauto S TG-N / Autosera S TG-N	99 *1	± 2. 1 *1
	Pureauto S TG-N / Autosera S TG-N	107. 3 *2	± 2. 0 *2
	Cholestest TG / QUALIGENT TG	107. 3 *2	± 2. 0 *2

\* 1) : Assigned value by using glycerol as the standard.

\* 2) : Assigned value by using the JSCC recommended method.

**SEKISUI MEDICAL CO.,LTD.**

1-3, Nihonbashi 2-chome, Chuo-ku, Tokyo 103-0027 JAPAN

"Cholestest", "QUALIGENT", "Autosera", and "Pureauto" are trademarks or registered trademarks of SEKISUI MEDICAL CO., LTD. in Japan and/or other countries.

All other brands, product names and service names are trademarks or registered trademarks of their respective companies.